

BCL2L14 antibody



Catalog Number: 108426

Product name

BCL2L14 antibody

Specificity

Human; other species not tested.

Antibody description

BCL2L14 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue. Positive WB detected in human lung tissue. Observed molecular weight by Western-blot: 37 kDa

Preparation

This antibody was obtained by immunization of BCL2L14 recombinant protein (Accession Number: NM_001370268). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:117099(BCL2L14 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:117099(BCL2L14 antibody) at dilution of 1:100 (under 10x lens)