

# BCMO1 antibody

Catalog Number: 108433

## Product name

BCMO1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

BCMO1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive WB detected in COLO 320 cells, mouse colon tissue. Observed molecular weight by Western-blot: 65 kDa

## Preparation

This antibody was obtained by immunization of BCMO1 recombinant protein (Accession Number: NM\_017429). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:117106(BCMO1 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:117106(BCMO1 Antibody) at dilution of 1:50 (under 40x lens)