

BCOR antibody

Catalog Number: 108435

Product name

BCOR antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

BCOR Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 186 kDa

Preparation

This antibody was obtained by immunization of BCOR recombinant protein (Accession Number: BC009675). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:117108(BCOR antibody) at dilution of 1:400