

BLVRA antibody

Catalog Number: 108455

Product name

BLVRA antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

BLVRA Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, L02 cells, mouse liver tissue. Positive IHC detected in human prostate cancer tissue, human pancreas tissue. Observed molecular weight by Western-blot: 37-40 kDa, 66 kDa

Preparation

This antibody was obtained by immunization of BLVRA recombinant protein (Accession Number: NM_000712). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117190(BLVRA antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:117190(BLVRA antibody) at dilution of 1:100 (under 10x lens)