

BMP10 antibody



Catalog Number: 108459

Product name

BMP10 antibody

Specificity

Human; other species not tested.

Antibody description

BMP10 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, MCF7 cells. Observed molecular weight by Western-blot: 48 kDa

Preparation

This antibody was obtained by immunization of BMP10 recombinant protein (Accession Number: NM_014482). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:117194(BMP10



antibody) at dilution of 1:600

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:117194(BMP10 antibody) at dilution of 1:600