

# BMP15 antibody



Catalog Number: 108460

## Product name

BMP15 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

BMP15 Rabbit Polyclonal antibody. Positive WB detected in A2780 cells, mouse ovary tissue. Positive IHC detected in human prostate cancer tissue. Observed molecular weight by Western-blot: 50 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_005448). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations




A2780 cells were subjected to SDS PAGE followed by western blot with Catalog No:117195(BMP15



antibody) at dilution of 1:300

mouse ovary tissue were subjected to SDS PAGE followed by western blot with Catalog No:117195(BMP15 antibody) at dilution of 1:300




Immunohistochemistry of paraffin-embedded human prostate cancer slide using Catalog No:117195(BMP15 Antibody) at dilution of 1:50 (under 10x lens)

# BMP15 antibody



Catalog Number: 108460

---



Immunohistochemistry of paraffin-embedded human prostate cancer slide using Catalog No:117195(BMP15 Antibody) at dilution of 1:50 (under 40x lens)