

# BMP4 antibody



Catalog Number: 108463

## Product name

BMP4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

BMP4 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, COLO 320 cells, human skeletal muscle tissue, rat heart tissue. Positive IHC detected in human lung tissue, human prostate hyperplasia tissue. Observed molecular weight by Western-blot: 50 kDa

## Preparation

This antibody was obtained by immunization of BMP4 recombinant protein (Accession Number: NM\_001347914). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


## Validations



Immunohistochemical of paraffin-embedded human lung using Catalog No:117198(BMP4 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human lung using Catalog No:117198(BMP4 antibody) at dilution of 1:50 (under 40x lens)




mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:117198(BMP4 Antibody) at dilution of 1:300

# BMP4 antibody



Catalog Number: 108463

---



rat heart tissue were subjected to SDS PAGE  
followed by western blot with Catalog  
No:117198(BMP4 Antibody) at dilution of 1:600