

BMP5 antibody



Catalog Number: 108464

Product name

BMP5 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

BMP5 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, HUVEC cells. Observed molecular weight by Western-blot: 52 kDa

Preparation

This antibody was obtained by immunization of BMP5 recombinant protein (Accession Number: NM_021073). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:117199(BMP5 Antibody) at dilution of 1:600