

BNC2 antibody



Catalog Number: 108472

Product name

BNC2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

BNC2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse uterus tissue. Observed molecular weight by Western-blot: 122-140 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_017637). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117207(BNC2 Antibody) at dilution of 1:600