

# BNIP1 antibody

Catalog Number: 108473

## Product name

BNIP1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

BNIP1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in human brain tissue, HeLa cells, human liver tissue, mouse skeletal muscle tissue, rat brain tissue. Observed molecular weight by Western-blot: 26 kDa

## Preparation

This antibody was obtained by immunization of BNIP1 recombinant protein (Accession Number: NM\_001205). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117208(BNIP1 antibody) at dilution of 1:500



IP Result of anti-BNIP1 (IP:Catalog No:117208, 3ug; Detection:Catalog No:117208 1:500) with mouse brain tissue lysate 6000ug.