

Beta galactosidase antibody



Catalog Number: 108482

Product name

Beta galactosidase antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Beta galactosidase Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells. Positive IF detected in Hela cells. Observed molecular weight by Western-blot: 67 kDa

Preparation

This antibody was obtained by immunization of Beta galactosidase recombinant protein (Accession Number: NM_000404). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:117127(GLB1 antibody) at dilution of 1:400



Immunofluorescent analysis of Hela cells, using GLB1 antibody Catalog No:117127 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).