

BHLHE41 antibody

Catalog Number: 108494

Product name

BHLHE41 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

BHLHE41 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human heart tissue, mouse testis tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 51 kDa, 75 kDa

Preparation

This antibody was obtained by immunization of BHLHE41 recombinant protein (Accession Number: NM_030762). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117139(BHLHE41 antibody) at dilution of



1:500
+ - - Normal Rabbit IgG
- + - BHLHE41 antibody

IP Result of anti-BHLHE41 (IP:Catalog No:117139, 4ug; Detection:Catalog No:117139 1:1000) with mouse brain tissue lysate 3440ug.