

BHMT2 antibody



Catalog Number: 108496

Product name

BHMT2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

BHMT2 Rabbit Polyclonal antibody. Positive WB detected in rat liver tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 40-45 kDa

Preparation

This antibody was obtained by immunization of BHMT2 recombinant protein (Accession Number: NM_017614). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



rat liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:117141(BHMT2 Antibody) at dilution of 1:1000