

# BIN3 antibody



Catalog Number: 108504

---

## Product name

BIN3 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

BIN3 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells, A375 cells, HeLa cells, mouse testis tissue. Observed molecular weight by Western-blot: 25kd

## Preparation

This antibody was obtained by immunization of BIN3 recombinant protein (Accession Number: XM\_047421995). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:117149(BIN3 antibody) at dilution of 1:500