

BIRC8 antibody



Catalog Number: 108505

Product name

BIRC8 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

BIRC8 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human brain tissue. Observed molecular weight by Western-blot: 30-35 kDa

Preparation

This antibody was obtained by immunization of BIRC8 recombinant protein (Accession Number: NM_033341). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117150(BIRC8 antibody) at dilution of 1:200



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117150(BIRC8 antibody) at dilution of 1:300