

ATXN3L antibody



Catalog Number: 108508

Product name

Polyclonal

ATXN3L antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

ATXN3L Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 41 kDa

Dilutions


Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of ATXN3L recombinant protein (Accession Number: NM_001135995). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:117069(ATXN3L Antibody) at dilution of 1:600