

AUP1 antibody

Catalog Number: 108511

Product name

AUP1 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

AUP1 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HEK-293 cells, Jurkat cells. Positive IP detected in HEK-293 cells. Positive IHC detected in human liver cancer tissue. Positive IF detected in A431 cells. Observed molecular weight by Western-blot: 40kd

Validations




Preparation

This antibody was obtained by immunization of AUP1 recombinant protein (Accession Number: BC033646). Purification method: Antigen affinity purified.

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:117072(AUP1 antibody) at dilution of 1:400

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:

WB: 1:500-1:5000


Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:117072(AUP1 antibody) at dilution of 1:100 (under 10x lens)



AUP1 antibody


Catalog Number: 108511

Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:117072(AUP1 antibody) at dilution of 1:100 (under 40x lens)



IP Result of anti-AUP1 (IP:Catalog No:117072,

3ug; Detection:Catalog No:117072 1:500) with HEK-293 cells lysate 880ug.



Immunofluorescent analysis of A431 cells using Catalog No:117072(AUP1 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG