

# AZI1 antibody



Catalog Number: 108521

## Product name

AZI1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

AZI1 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 131 kDa

## Preparation

This antibody was obtained by immunization of AZI1 recombinant protein (Accession Number: XM\_047435657). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:117082(AZI1 Antibody) at dilution of 1:1000



Immunofluorescent analysis of HeLa cells using Catalog No:117082(AZI1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)