

B3GALNT1 antibody



Catalog Number: 108524

Product name

B3GALNT1 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

B3GALNT1 Rabbit Polyclonal antibody. Positive WB detected in mouse placenta tissue, human placenta tissue, mouse brain tissue, mouse colon tissue. Positive IHC detected in human heart. Observed molecular weight by Western-blot: 40-45kd

Preparation

This antibody was obtained by immunization of B3GALNT1 recombinant protein (Accession Number: NM_001349138). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

Validations




mouse placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:117085(B3GALNT1 antibody) at dilution of



1:1000

Immunohistochemical of paraffin-embedded human heart using Catalog No:117085(B3GALNT1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:117085(B3GALNT1 antibody) at dilution of 1:50 (under 40x lens)