

B3GALNT2 antibody



Catalog Number: 108525

Product name

B3GALNT2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

B3GALNT2 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Positive IHC detected in human heart. Observed molecular weight by Western-blot: 55-60kd

Preparation

This antibody was obtained by immunization of B3GALNT2 recombinant protein (Accession Number: NM_152490). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:117086(B3GALNT2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human heart using Catalog No:117086(B3GALNT2 antibody) at dilution of 1:100 (under 40x lens)