

# B4GALNT1 antibody



Catalog Number: 108533

## Product name

B4GALNT1 antibody

## Specificity

Human; other species not tested.

## Antibody description

B4GALNT1 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, A375 cells. Observed molecular weight by Western-blot: 67.5kd

## Preparation

This antibody was obtained by immunization of B4GALNT1 recombinant protein (Accession Number: NM\_001276469). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:117094(B4GALNT1 antibody) at dilution of 1:400