

BTBD6 antibody

Catalog Number: 108537

Product name

BTBD6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

BTBD6 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Positive IP detected in Y79 cells. Observed molecular weight by Western-blot: 53-58 kDa

Preparation

This antibody was obtained by immunization of BTBD6 recombinant protein (Accession Number: BC042525). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:108537(BTBD6 Antibody) at dilution of 1:300



IP Result of anti-BTBD6 (IP:Catalog No:108537, 4ug; Detection:Catalog No:108537 1:300) with Y79 cells lysate 800ug.