

BTK antibody

Catalog Number: 108544

Product name

BTK antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

BTK Rabbit Polyclonal antibody. Positive IP detected in mouse spleen tissue. Positive WB detected in mouse spleen tissue, L02 cells. Observed molecular weight by Western-blot: 70-76 kDa

Preparation

This antibody was obtained by immunization of BTK recombinant protein (Accession Number: NM_000061). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:108544(BTK antibody) at dilution of 1:300



IP Result of anti-BTK (IP:Catalog No:108544, 4ug; Detection:Catalog No:108544 1:300) with mouse spleen tissue lysate 4800ug.