

BTN2A1 antibody



Catalog Number: 108545

Product name

BTN2A1 antibody

Specificity

Human; other species not tested.

Antibody description

BTN2A1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, MCF7 cells, mouse spleen tissue, rat spleen tissue. Observed molecular weight by Western-blot: 69kd

Preparation

This antibody was obtained by immunization of BTN2A1 recombinant protein (Accession Number: NM_078476). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108545(BTN2A1 antibody) at dilution of 1:600