

C10orf125 antibody



Catalog Number: 108561

Product name

C10orf125 antibody

Specificity

Human; other species not tested.

Antibody description

C10orf125 Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer tissue.

Preparation

This antibody was obtained by immunization of C10orf125 recombinant protein (Accession Number: NM_198472). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:108561(C10orf125 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:108561(C10orf125 Antibody) at dilution of 1:50 (under 40x lens)