

C10orf58 antibody



Catalog Number: 108565

Product name

C10orf58 antibody

Immunogen

[Human C10orf58 Recombinant protein \(GST tag\)](#)

Specificity

Human, Rat; other species not tested.

Antibody description

C10orf58 Rabbit Polyclonal antibody. Positive WB detected in ROS1728 cells. Positive IF detected in ROS1728 cells. Observed molecular weight by Western-blot: 24-30 kDa

Preparation

This antibody was obtained by immunization of C10orf58 recombinant protein (Accession Number: NM_001243779). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

Validations



ROS1728 cells were subjected to SDS PAGE followed by western blot with Catalog No:108565(C10orf58 Antibody) at dilution of



1:300

Immunofluorescent analysis of ROS1728 cells using Catalog No:108565(C10orf58 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)