

Brachyury-T antibody



Catalog Number: 108574

Product name

Brachyury-T antibody

Specificity

human, mouse; other species not tested.

Antibody description

Brachyury-T Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 44 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003181). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:117221(Brachyury-T antibody) at dilution of 1:500