

BRCA2 antibody



Catalog Number: 108578

Product name

BRCA2 antibody

Specificity

Human, Rat; other species not tested.

Antibody description

BRCA2 Rabbit Polyclonal antibody. Positive IHC detected in rat brain tissue. Positive WB detected in human brain tissue, HeLa cells. Observed molecular weight by Western-blot: 384kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000059).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117225(BRCA2 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded rat brain using Catalog No:117225(BRCA2 antibody) at dilution of 1:50 (under 10x lens)