

BRD8 antibody

Catalog Number: 108583

Product name

BRD8 antibody

Specificity

Human; other species not tested.

Antibody description

BRD8 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, human brain tissue. Observed molecular weight by Western-blot: 120 kDa

Preparation

This antibody was obtained by immunization of BRD8 recombinant protein (Accession Number: BC008039). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117230(BRD8 antibody) at dilution of 1:500



IP Result of anti-BRD8 (IP:Catalog No:117230, 4ug; Detection:Catalog No:117230 1:300) with HeLa cells lysate 1200ug.