

Brevican antibody



Catalog Number: 108586

Product name

Brevican antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Brevican Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, COLO 320 cells, SH-SY5Y cells, Y79 cells. Observed molecular weight by Western-blot: 145 kDa, 80 kDa

Preparation

This antibody was obtained by immunization of Brevican recombinant protein (Accession Number: BC010571). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117233(BCAN antibody) at dilution of 1:500