

BRIP1 antibody



Catalog Number: 108588

Product name

BRIP1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

BRIP1 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 100 kDa

Preparation

This antibody was obtained by immunization of BRIP1 recombinant protein (Accession Number: NM_032043). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:117235(BRIP1 Antibody) at dilution of 1:600