

BRUNOL5 antibody



Catalog Number: 108599

Product name

BRUNOL5 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

BRUNOL5 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas. Positive WB detected in mouse skeletal muscle tissue, mouse testis tissue. Observed molecular weight by Western-blot: 52kd

Preparation

This antibody was obtained by immunization of BRUNOL5 recombinant protein (Accession Number: XM_006722832). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:117246(BRUNOL5 antibody) at dilution of



1:1500
Immunohistochemical of paraffin-embedded human gliomas using Catalog No:117246(BRUNOL5 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:117246(BRUNOL5 antibody) at dilution of 1:50

BRUNOL5 antibody



Catalog Number: 108599

(under 40x lens)