

# C16orf61 antibody



Catalog Number: 108606

## Product name

C16orf61 antibody

## Immunogen

[Human C16orf61 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

C16orf61 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, human kidney tissue. Observed molecular weight by Western-blot: 9 kDa

## Preparation

This antibody was obtained by immunization of C16orf61 recombinant protein (Accession Number: NM\_020188). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108606(C16orf61 Antibody) at dilution of 1:600