

# C17orf85 antibody



Catalog Number: 108613

## Product name

C17orf85 antibody

## Specificity

Human; other species not tested.

## Antibody description

C17orf85 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A549 cells, COLO 320 cells, HEK-293 cells, Jurkat cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 71 kda

## Preparation

This antibody was obtained by immunization of C17orf85 recombinant protein (Accession Number: NM\_018553). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108613(C17orf85 Antibody) at dilution of 1:1000



Immunofluorescent analysis of HeLa cells using Catalog No:108613(C17orf85 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG