

C19orf21 antibody



Catalog Number: 108619

Product name

C19orf21 antibody

Specificity

Human; other species not tested.

Antibody description

C19orf21 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, HepG2 cells. Observed molecular weight by Western-blot: 75 kDa

Preparation

This antibody was obtained by immunization of C19orf21 recombinant protein (Accession Number: XM_011527685). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:108619(C19orf21 Antibody) at dilution of 1:2000