

# C1orf116 antibody



Catalog Number: 108633

---

## Product name

C1orf116 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

C1orf116 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, mouse testis tissue. Observed molecular weight by Western-blot: 64 kDa, 38 kDa

## Preparation

This antibody was obtained by immunization of C1orf116 recombinant protein (Accession Number: NM\_023938). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108633(C1orf116 antibody) at dilution of 1:800