

# C11orf54 antibody



Catalog Number: 108640

## Product name

C11orf54 antibody

## Specificity

Human, Rat; other species not tested.

## Antibody description

C11orf54 Rabbit Polyclonal antibody. Positive IHC detected in mouse liver tissue. Positive WB detected in rat liver tissue. Observed molecular weight by Western-blot: 31 kDa

## Preparation

This antibody was obtained by immunization of C11orf54 recombinant protein (Accession Number: BC007110). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



rat liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:108640(C11orf54 Antibody) at dilution of



1:300  
Immunohistochemistry of paraffin-embedded mouse liver tissue slide using Catalog No:108640(C11orf54 Antibody) at dilution of 1:50 (under 40x lens)