

C11orf79 antibody



Catalog Number: 108645

Product name

C11orf79 antibody

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

C11orf79 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human skin cancer tissue. Positive WB detected in A549 cells, HeLa cells, PC-3 cells. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of C11orf79 recombinant protein (Accession Number: NM_017841). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:1000


Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:108645(C11orf79 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:108645(C11orf79 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:108645(C11orf79 Antibody) at dilution of 1:50

C11orf79 antibody



Catalog Number: 108645

(under 10x lens)