

C13orf1 antibody



Catalog Number: 108651

Product name

Polyclonal

C13orf1 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

C13orf1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 17 kDa, 20 kDa

Dilutions


Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of C13orf1 recombinant protein (Accession Number: NM_001127482). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108651(C13orf1 Antibody) at dilution of 1:600

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality