

# C13orf33 antibody



Catalog Number: 108654

## Product name

C13orf33 antibody

## Specificity

Human; other species not tested.

## Antibody description

C13orf33 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Positive IP detected in HepG2 cells. Observed molecular weight by Western-blot: 34-45 kDa

## Preparation

This antibody was obtained by immunization of C13orf33 recombinant protein (Accession Number: NM\_032849). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:108654(C13orf33 Antibody) at dilution of 1:600

IP Result of anti-C13orf33 (IP:Catalog No:108654, 4ug; Detection:Catalog No:108654 1:300) with HepG2 cells lysate 4000ug.