

C14orf156 antibody



Catalog Number: 108659

Product name

C14orf156 antibody

Specificity

Human; other species not tested.

Antibody description

C14orf156 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HeLa cells, MDA-MB-453s cells. Observed molecular weight by Western-blot: 12 kDa

Preparation

This antibody was obtained by immunization of C14orf156 recombinant protein (Accession Number: NM_031210). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108659(C14orf156 Antibody) at dilution of 1:600