

C14orf159 antibody

Catalog Number: 108660

Product name

C14orf159 antibody

Immunogen

[Human C14orf159 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

C14orf159 Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue.

Preparation

This antibody was obtained by immunization of C14orf159 recombinant protein (Accession Number: NM_001102367). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations

Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:108660(C14orf159 Antibody) at dilution of 1:50 (under 40x lens)