

# C1orf163 antibody



Catalog Number: 108672

## Product name

C1orf163 antibody

## Immunogen

[Human C1orf163 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

C1orf163 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 25-28 kDa

## Preparation

This antibody was obtained by immunization of C1orf163 recombinant protein (Accession Number: NM\_023077). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108672(C1orf163 Antibody) at dilution of 1:1000