

# C1orf26 antibody

Catalog Number: 108675

## Product name

C1orf26 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

C1orf26 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 100-105 kDa

## Preparation

This antibody was obtained by immunization of C1orf26 recombinant protein (Accession Number: NM\_017673). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:200-1:2000

## Validations



IP Result of anti-C1orf26 (IP:Catalog No:108675, 5ug; Detection:Catalog No:108675 1:500) with mouse testis tissue lysate 8000ug.