

C1orf27 antibody



Catalog Number: 108676

Product name

C1orf27 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

C1orf27 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of C1orf27 recombinant protein (Accession Number: XM_011509660). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108676(C1orf27 Antibody) at dilution of 1:300