

# C1QA antibody

Catalog Number: 108685

## Product name

C1QA antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

C1QA Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in human spleen tissue, human colon tissue, human heart tissue, human liver tissue, human lung tissue, mouse liver tissue, mouse serum tissue. Observed molecular weight by Western-blot: 30 kDa

## Preparation

This antibody was obtained by immunization of C1QA recombinant protein (Accession Number: NM\_015991). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:108685(C1QA antibody) at dilution of 1:500



IP Result of anti-C1QA (IP:Catalog No:108685, 4ug; Detection:Catalog No:108685 1:300) with mouse liver tissue lysate 4000ug.