

C1QTNF3 antibody

Catalog Number: 108688

Product name

C1QTNF3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

C1QTNF3 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human small intestine tissue.

Preparation

This antibody was obtained by immunization of C1QTNF3 recombinant protein (Accession Number: BC016021). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:108688(C1QTNF3 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:108688(C1QTNF3 Antibody) at dilution of 1:50 (under 40x lens)