

C1R antibody



Catalog Number: 108689

Product name

Polyclonal

C1R antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

C1R Rabbit Polyclonal antibody. Positive WB detected in human blood tissue. Observed molecular weight by Western-blot: 85 kDa, 30 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of C1R recombinant protein (Accession Number: NM_001733). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:108689(C1R Antibody) at dilution of 1:1000