

C20orf152 antibody



Catalog Number: 108694

Product name

C20orf152 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

C20orf152 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue. Observed molecular weight by Western-blot: 67-70 kDa

Preparation

This antibody was obtained by immunization of C20orf152 recombinant protein (Accession Number: NM_001207076). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:108694(C20orf152 Antibody) at dilution of 1:300