

# C20orf194 antibody



Catalog Number: 108697

## Product name

C20orf194 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

C20orf194 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, rat testis tissue. Observed molecular weight by Western-blot: 132 kd

## Preparation

This antibody was obtained by immunization of C20orf194 recombinant protein (Accession Number: BC122529). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108697(C20orf194 Antibody) at dilution of 1:1000